

Head lice – a growing problem

Head lice have become a growing problem during the last ten years. It is mostly children between the ages of 3-11 years that are affected. Usually the problem is noted at the beginning of each school term, but infestations occur all year long. Head lice can neither fly nor jump. To get lice it is necessary to have close hair contact with an infested person. Because people today have more social contacts the risk of getting infested with head lice has increased. Head lice is a growing problem not only among children but among whole families, since lice easily spread to all family members before the infestation is discovered and treated.

Why should I use the NitFree lice comb?

By using the NitFree lice comb, you can remove both lice and nits at the same time, disrupting the lice's biological lifecycle. Nits are not able to hatch into new lice. When removing nits – no new lice will be hatching.

Läkemedelsverket highlights the importance of using a lice comb in all head lice treatments. What distinguishes NitFree from other products is NitFree's patented micro grooved teeth, which enables it to simultaneously remove both lice and nits. Nits are difficult to remove from the hair, as these are firmly attached to the hair with a glue like subject. The distance between NitFree's teeth is less than 0,2 mm, making it possible to trap, destroy and remove even the smallest nits while combing.

NitFree is made of stainless steel. It is durable and both appearance and shape is kept despite repeated use, since NitFree is made with precision based laser technology. The comb can be sterilized and used by all family members without risk of further infestation.

NitFree is particularly suitable for sensitive people or persons with allergies, is environmentally friendly, does not contribute to resistance to insecticides among lice, and can be used proactively to you quickly discover and take control of head lice.

To really succeed in getting rid of lice and continue to be lice-free, continue using a NitFree. Combing is the most important part of treating lice no matter what treatment you choose. If you want to use pediculicides for killing lice, you can also use them together with NitFree

What makes NitFree unique?

Durable and possible to sterilize

Made of stainless steel, maintains both shape and form.

Non-slip grips

Easy to use, won't slip with wet or soapy hands.

Long teeth

Gives better reach. More nits and lice can be removed with each comb stroke.

Exact distance between comb teeth


Teeth are separated from each other by a carefully-calculated gap, smaller than even the tiniest nit. Even the very smallest nits are removed.



Unique micro grooved teeth

Nits are destroyed and removed without pulling, cutting or hurting the hair they cling to.

Rounded ends, no sharp points to prick, scratch or bruise the scalp.



NitFree testimonials

"We struggled for months to get rid of all lice and eggs, with no success. When I read what praise NitFree got (I am in marketing and sales) I thought, "well let's see if it works.." Well after using it – IT WORKS! We've combed and combed, and now we've removed those eggs that we earlier had to remove with tweezers, only by a simple stroke. If you doubt whether or not to buy NitFree here is my strong recommendation to BUY!"

Sofie, mother of two "thick-haired" children

"Hi, I wanted to thank you for your fantastic comb! After less than 2 weeks we are now free from lice. First I doubted it would work on our afro style hair, but it worked really well. I even think NitFree is easier to use than the combs we usually use for combing our afros!"

Tanja, mother of two - all with afrostyle hair

"We got rid of all lice in only 3 days, without using any chemical lice killing product. We are very pleased with NitFree. It's great to always have an easy solution for the problem right at home."

Per, father of three children

How to use NitFree lice comb

Use NitFree in the bathroom, then you can easily clean the comb under running water between each comb stroke.

1. Keep the hair dry or wet it lightly to avoid static electricity.

2.  Use a regular comb or brush to untangle the hair.

3.  Divide the hair into smaller sections and use NitFree to carefully comb through each section. It is important to comb as close to the scalp as possible to remove all nits. To better see what is removed you can strike the comb lightly on a white piece of paper or a mirror. The sections that have been combed through should be kept separate from those not yet combed.

4.  Clean NitFree between each comb stroke. Carefully pull the teeth slightly apart under running water to remove all lice and nits.

5. Sterilize NitFree when you have finished combing by putting it in boiling water for 1 minute. Don't forget to clean any other combs or brushes you might have used to untangle the hair.

Repeat the procedure above on a daily basis until you get rid of all lice and nits. We also recommend that you check if other members of your family have lice or nits to avoid further contamination.

We recommend that you continue to use NitFree regularly (once a week) to detect problems early and avoid the further spread of lice.

www.LusFri.nu

NitFree is available in healthcare shops and at: www.LusFri.nu.

For more information and instructions for use in different languages, please enter: www.LusFri.nu or send an email to: info@lusfri.nu.

Sales and distribution in the Nordic countries:

LusFri

Tutviksvägen 15
136 75 Haninge
Tel: 073 - 503 3436

www.LusFri.nu

Head lice?

Get LusFri.nu

Take control – eliminate lice quickly by starting with the nits.



NitFree – lice comb with unique micro grooved teeth which removes lice and eggs at the same time.

Strong and durable. Recommended for allergic and sensitive persons.